

# Industrial Evaporation

## The next generation



# Quality that has impressed for decades

For more than 50 years the rotary evaporators from Büchi Labortechnik AG have focused on ease of use, durability and quality. For this reason, 20-year-old instruments can still be found in pharmaceutical and chemical laboratories all over the world. Such expectations to meet high quality standards are also found within industrial rotary evaporation line. For these reasons and Büchi's promise to continue to deliver quality instrumentation, trust that the processing of your valuable products will deliver unbeatable results with the Büchi Rotavapor.





# The perfect Rotavapor for every laboratory

3

The industrial Rotavapors were developed to be of universal application, covering an almost endless range of applications. Most typically the industrial Rotavapors are used in scale-up processes.

## Typical applications for the industrial Rotavapors are:

- Concentrations in kilo scale as well as in production
- Recrystallizations
- Drying
- Solvent work-up (cleaning)
- Extractions of plant products

### R-220 SE

The Rotavapor R-220 SE (Second Edition) is a further development of the industry approved Rotavapor R-220. Many proven elements have been carried over from the original R-220. However, a few new features have been introduced in relation to safety and ease of use.



### R-220 EX

The ATEX-compliant version of the R-220 was specially developed for applications in EX zones and is compliant with the extremely stringent ATEX 95 directives. The Rotavapor R-220 EX is suitable for EX zones class 1 and 2, and provides the highest possible safety combined with maximum convenience.



### R-250

Up to 30 liters of solvent can be quickly and safely distilled in the 50-liter evaporating flask for the R-250, and that at a performance of up to 30 liters/hour (applies to acetone).



### R-250 EX

For the distillation of large volumes in EX-protected pharmaceutical production facilities, the R-250 EX is the optimum tool. For this environment, FDA-compliant materials are a necessity.



# Easy, safe, automated processes

The Buchi Rotavapors feature intuitive operation. Numerous technical features make the industrial Rotavapors safe, robust instruments for continuous use in very demanding environments.



## Title Feed valve

Quick or continuous addition of sample without breaking vacuum



## Patented EasyClamp connections

Secures glass joints quickly, easily and safely



## Patented flask clamp

Snap action mechanism – quick, easy and safe



## Spring loaded shut-off valves

Prevents glass breakage



## 6" LCD display

Large informative digital display at a glance



## Manual flask handler

The safest way to carry and mount flasks




## Integrated vacuum controller

Easy and reproducible process control




## Universal heating bath

For use with water or oil



## Coated glass

All units equipped with P+G glass providing maximum user protection



## EX protection

EX versions comply with ATEX 95 directives



### Warning alarm

Floating connection for external alarm equipment such as warning buzzer

**NEW**



### Level sensors

No more overfilling of the receiving flask

### Foam detector

Unsupervised operation with foaming products



### Integrated 230 V power supply

to connect the Vacuum pump V-710

**NEW**



### Data storage

Easy data storage on USB stick



**NEW**

**NEW**



### FDA compliant materials

Maximum product assurance

**NEW**



### Protective enclosure

Additional safety for the instrument

### Custom units

Customer specific solutions from Buchi



## Other impressive characteristics:

### Safe and easy operation

- Direct access to all glass components
- PTFE base outlet valve for chemical resistant easy emptying
- Front operated bath replenishment to prevent bath from running dry
- Graduated receiving flasks for measurable reproducible results
- Soft starting of the rotation for safer operation

### Process optimization

- Distillate condenser on all reflux glass assemblies to cool the distillate and reduce to loss of solvent.
- Optional vapor duct with integrated frit prevents glass assembly contamination during powder drying

### Future innovation

- RS485 and USB connection for the integration of future accessories

You will find detailed information via Internet at [www.buchi.com](http://www.buchi.com).

# Increase efficiency

**Our industrial rotary evaporators were developed with a focus on efficiency and safety.**

**To increase efficiency there are generally two options. Either increase the distillation capacity through more heating capacity or optimize operation and monitoring to minimize personnel costs. Buchi rotary evaporators have considered and implemented both approaches to increase process efficiency.**

## **Distillation capacity**

The optimal combination of heating capacity and cooling area allows for distillation capacities up to 30 liters of acetone/hour on the R-250 (EX). The R-220 SE also achieves a remarkable 19 liters/hour.

## **Operation and monitoring**

The following features make working with the large rotary evaporators easier such that the efficiency and therefore the cost savings is increased.



### **Patented flask attachment with snap action mechanism together with manual flask handler**

The evaporating flask can be fastened easily and safely to the Rotavapor by only one person. As a fitting is not used, the following advantages are obtained:

- Quick: – No additional accessory needed to secure joint connection.
- Safe: – No self loosening of the fitting.
  - No fastening parts that can be lost, as the flask attachment is part of the Rotavapor.
  - Safe handling of the flask during transport and during fastening, as the flask cannot roll away.



### **Patented EasyClamp connections on glass joints**

The quickest and, at the same time, the safest way to attach and remove glass components.

- Safe: – Very robust material and gentle on glass.
- Quick: – Opening and closing with only 1 screw.



### Vacuum controller

The integrated vacuum controller on the R-220 SE regulates and optimizes the distillation process with minimal user involvement and supervision. Functions such as Timer, Gradient and Repeat are only three of the special features.

- Safe: – Stable conditions and reproducible processes.
- Dependable: – Tightly controlled settings result in reproducible processes.

**NEW**

## Optional accessories for time and cost savings



### Foam detector

The internal PTFE foam sensor detects and breaks down rising foam with a short aeration pulse. Therefore, you can use the foam sensor to eliminate the time consuming monitoring process of extract reduction.

- Efficient: – Reduction in monitoring results in cost savings.



### Level sensors

Adjustable receiving flask level sensors can monitor various fill levels. As soon as they are full, the Rotavapor switches to a safe state (bath lowered, heater and rotation off).

- Safe: – No overfilling of the receiving flask.
- Dependable: – Measurable concentration.
- Efficient: – Reduced monitoring effort.

**NEW**



### Data storage

Particularly in the pharmaceutical industry, saving process data is very important. The USB module is an option available to save all relevant process data to a USB stick. This data can then be transferred and saved to a PC.

- Safe: – Subsequent process evaluation.
- Efficient: – Process monitoring without personnel costs.

**NEW**



# Safety without compromises

**There are various features that make the Rotavapor a safe tool for daily distillation tasks in both the laboratory and production.**



## Product assurance

Maximum protection for thermosensitive products against excessively high temperatures.

- An adjustable upper temperature limit prevents unintentional change of the setpoint.
- An overtemperature sensor switches the Rotavapor off when setpoint temperature has been exceeded by 15 °C.
- A mechanical over temperature protection feature protects the heater from burning out.

The EX versions have level monitoring which prevent the baths from running dry.

## Automatic bath lowering

In the event of a problem, the bath is lowered automatically, the heater switches off and the rotation stops.



## Level monitoring for the receiving flasks (optional)

If the receiving flask is full as indicated by the level sensor, the R-220 SE switches off.

**NEW**

## Process reliability

Vacuum controller in combination with speed controlled vacuum pump, coolant monitoring and overall precise regulation of the evaporation parameters guarantee reproducible results.



## User protection

Protection of the operator is of high value with all Buchi instrumentation. This means that even the simplest and most economical Rotavapor offers maximum protection, without having to purchase additional safety equipment.

Standard features include:

- “Plastic+Glas” coating on all glass components with the exception of evaporating flasks and EX versions. This coating offers maximum protection for personnel against implosions and glass breakage without hindering operation.
- Spring-loaded shut-off valves prevent glass breakage as a result of excessive force during closing.
- EasyClamp connections for quick and safe assembly of glass joints .
- Patented evaporation flask clamp allows safe single user operation.



## Instrument protection


Additional protection of the instruments is often required in production facilities. For this reason a safety shield made of PMMA is available for the R-220 SE; this shield protects the entire glass assembly.

- Unrestricted visibility of the process.
- Front access to the controls via doors.
- The coated glass components ensure the protection of the user even with the doors open.

The pivoting polycarbonate bath shield protects, on the one hand, the evaporating flask and, on the other hand, the user against splashing.

- Transparent shield for continuous visibility of the product throughout the entire evaporation process.
- Pneumatic spring for easy pivoting and retention of the shield.
- Safety feature – rotation stops when shield is open.
- Optional stainless steel cover protects the evaporation flask from above.


A steel tube frame is also available as an accessory that provides additional side protection for the R-220 SE.



# The Rotavapor® R-220 SE is setting new standards

The Rotavapor R-220 SE (Second Edition) is based on the extremely successful R-220.


This new edition R-220 improves upon safety standards and flexibility.



**The integrated vacuum controller pre-programmed with numerous functions, making distillations quicker, more reproducible and safer.**

**The most important features include:**

- Large display for clear indication of important data
- Control knob allows for simple menu operation
- Accurate pressure control
- Prepare and save gradients for reproducible distillations
- Timer function stops rotation, bath temperature and lift position after a programmed time is achieved
- Solvent library incorporates bath and coolant temperatures for optimal distillation
- Hysteresis free vacuum settings reduce the occurrence of bumping and foaming




**NEW**

**The large, clearly laid out display makes it possible to monitor the process at a glance. All parameters and settings are displayed to eliminate submenu navigation. These features include:**

- Actual bath temperature
- Vapor temperature
- Coolant temperature (requires optional sensor)
- Cooling status ON/OFF (requires optional sensor)
- Rotation speed (rpm)
- Heating status (ON/OFF)
- Bath shield (OPEN/CLOSED)
- Error messages
- Operating states in the continuous version (filling, emptying, automatic)
- Receiving flask level sensor status
- USB stick data storage activation (RECORD)

# EX protection in accordance with ATEX 95

Büchi's R-220/250 EX Rotavapors comply with the international standards (IEC and CENELEC) as well as the EC directive 94/9 (ATEX 95). The R-220 EX and R-250 EX Rotavapors support the pharmaceutical and chemical industries, where production under EX protected conditions is often the standard.



## Special features of the EX versions

- Hydraulic lift and bath rotation.
- Compliance with the EC directive 94/9 (ATEX 95).
- Maximum bath temperature is 150°C.
- Includes EX protected vacuum controller.
- All glass components are uncoated due to possible static charge.

## Zone classification (workplaces as per ATEX 137)

Europe: zones 1 and 2 (category 2 G), gas group IIC,  
Temperature class T3 or T4 (corresponds in the USA to:  
division 1, class 1, group C and D, temperature class T3/T4).

## User friendly operation also on the EX instruments



R-220 EX



R-250 EX with vacuum controller




# The continuous option – for uninterrupted operation

The continuous option allows for the distillation or concentration of several hundred liters of product when using the Rotavapor R-220 SE, R-250 or R-250 EX. These instruments are particularly suitable for solvent recycling or concentrations in industrial preparative chromatography.

The principle could not be easier and does not require any programming. The control system takes over the following functions:

- Automatically empties the receiving flask when full.
- Simultaneously refills the evaporation flask.
- Monitors fill levels in both the evaporating and receiving flasks.
- Controls PTFE fill and drain valves.
- Safety shut down in case of flooding the system, empty product tank or blocked drain.

- **Quick**  
Continuous concentration in preparative chromatography.
- **Automatic**  
Solvent recycling for large volumes.
- **Flexible**  
Manual and/or continuous operation is possible
- **Safe**  
Monitored with sensors.
- **Dependable**  
Easy operation without programming.




# Maximum performance with the right glass assembly


Numerous glass assemblies and accessories provide a solution for almost every application. Optimal combinations are provided below.

R-220 SE/EX	D	D2	DB2	DB	RB	R	C
R-250/EX		D2	D3			RB2	R2
Use	Low boiling points and/or foaming products				High boiling points		Very low boiling point
		Minimum emissions			Reflux reactions		
			Reduced height				

Glass assemblies for the R-220 line




Glass assemblies for the R-250 line




Application examples

Application	Example	Glass assemblies	Accessory
Purifying contaminated ethanol	Hospital	D	
Concentrating foaming product	Plant extract	D	Foam sensor
Thermo sensitive product	Concentration	D2 at low temperature (possibly glass assy. C)	Coolant temperatur sensor
Concentration of large volumes	Industrial chromatography	Continuous version	
Recrystallization		R	
Drying		R	Vapor duct with frit

Evaporating performance



- 20-liter flask, acetone, 130 rpm
- 10-liter flask, acetone, 130 rpm
- 20-liter flask, acetone, 30 rpm



- 50-liter flask, acetone, 100 rpm
- 20-liter flask, acetone, 100 rpm
- 50-liter flask, acetone, 130 rpm

For other solvents the following maximum general figures apply:

Solvent	Water	Methanol	Ethanol	Toluene	Ethylacetate	Hexane	Trichlorethane
Evaporating performance (liter per hour) R-220 SE	4.0	9.5	13.5	25.0	26.0	38.0	28.0
Evaporating performance (liter per hour) R-250	5.8	14.0	19.5	36.5	38.0	55.0	40.5

Temperature difference between bath temperature and boiling point (vapor temperature) = 40 °C

# Accessories R-220 SE and R-220 EX

## Accessories for the R-220 SE



### Vacuum pump V-710

Speed controlled 3-stage vacuum pump with 4 diaphragm heads, high flow rate of 3.1 m³/h and low ultimate vacuum of 2 mbar.

	Order no.
Vacuum pump V-710	72000
with Woulff bottle and secondary condenser	72011



### Trolley

Trolley made of stainless steel for easy transportation of Rotavapor with shelf space for V-710 vacuum pump. Dimensions W x H x D 850 x 560 x 420 mm.

	Order no.
Trolley for R-220 SE and EX	41257



### Manual flask handler

Assists in easy mounting and removal of large volume flask along with safe transport and sturdy base for positioning.

	Order no.
Manual flask handler for 20-liter flask for R-220 SE and EX	41400



### Evaporating and drying flasks

Large flange (inner diameter 120 mm) allows for easy sample removal and cleaning.

	Order no.
6-liter flask	27470
10-liter flask	27469
20-liter flask	27468
10-liter drying flask	28592
20-liter drying flask	28593



### Foam detector

Internal sensor detects rising foam and triggers a short aeration pulse, eliminating foam

	Order no.
Foam sensor with distribution head cpl.	11056083



### Vacuum valve

Electrical valve for vacuum regulation when working with house vacuum system or non-Buchi vacuum pump.

	Order no.
Vacuum valve with 10 mm connections	11055928



### Vapor duct

A special vapor duct with frit P3 protects the condenser assembly against powder and dust during the drying process.

	Order no.
Vapor duct with frit P3	41100



### Flow monitor

Checks the flow of coolant, stopping operation when flow of coolant is insufficient or interrupted.

	Order no.
Flow monitor	11055971



### Coolant temperature sensor

Displays coolant temperature for optimal monitoring of distillation parameters.

	Order no.
Coolant temperature sensor	11055988



### Cooling water valve

Eliminates unnecessary water waste by stopping cooling water flow when not in use.

	Order no.
Cooling water valve	41191



### Level sensor

Receiving flask level sensor terminates operation when predetermined waste levels are achieved. Particularly helpful for concentration of product.

	Order no.
Level sensor	11056192



### USB module

Saves parameters from application runs (temperatures, pressure, rotational speed) on a USB stick for data transfer.

	Order no.
USB module cpl.	11056175



Safety shield for the following glass assemblies, with *one* receiving flask:  
D, D2, R, C

	Order no.
Standard safety shield	11055796

Safety shield for the following glass assemblies, with *two* receiving flask:  
D, D2, R, C

	Order no.
standard safety shield, two flasks	11055797

Safety shield for the following glass assemblies, with *one* receiving flask:  
DB, DB2, RB

	Order no.
Safety shield, Bullfrog	11055798

Safety shield for the following glass assemblies, with *two* receiving flask:  
DB, DB2, RB

	Order no.
Safety shield, Bullfrog, two flasks	11055799



Separate bath shield complete

	Order no.
Bath shield cpl.	11055364

Protects the glass assembly and the bath against falling objects (only in conjunction with corresponding safety shield)

	Order no.
Protective grid cpl.	11056081

Reinforced steel tube frame

	Order no.
Bullbar cpl.	11056082

### Accessories for the EX protected R-220 EX



#### ATEX compliant diaphragm vacuum pump

EX area compliant pump achieves a flow rate of 3.7 m<sup>3</sup>/h and ultimate vacuum of 3 mbar. Complete with separator (not sure what this is) and secondary condenser.

Order no.	
Vacuum pump MD 4C EX 230 V/50 Hz	40197



#### Trolley

Trolley made of stainless steel for easy transportation of Rotavapor with shelf space for V-710 vacuum pump. Dimensions W x H x D 850 x 560 x 420 mm.

Order no.	
Trolley for R-220 SE and EX	41257



#### Manual flask handler

Assists in easy mounting and removal of large volume flask along with safe transport and sturdy base for positioning.

Order no.	
Manual flask handler for 20-liter flask for R-220 SE and EX	41400



#### Evaporating and drying flasks

Large flange (inner diameter 120 mm) allows for easy sample removal and cleaning.

	Order no.
6-liter flask	27470
10-liter flask	27469
20-liter flask	27468
10-liter drying flask	28592
20-liter drying flask	28593



#### Safety shield

Pivoting stainless steel shield protects both user and flask

Order no.	
Safety shield for the EX area	46431



**Only accessories with this symbol are to be used in EX zones.**

# Accessories R-250 and R-250 EX

## Accessories for the Rotavapor R-250/R-250 EX



### Vacuum controller V-850/855

Display, control and regulation of vacuum. Includes solvent library and timer function, with additional gradient mode on the V-855.



### Vacuum valve

Electrical valve for vacuum regulation when working with house vacuum system or non-Buchi vacuum pump.



### Vacuum pump Sogevac SV 40

Rotary vane oil vacuum pump with gas ballast valve. Flow rate 40 m<sup>3</sup>/h with an ultimate pressure of < 2 mbar.



### Foam detector

Internal sensor detects rising foam and triggers a short aeration pulse, eliminating foam

Order no.	
Vacuum controller V-850 for the R-250	47293
Vacuum controller V-855 for the R-250	47292

Order no.	
Vacuum valve for R-250	31355

Order no.	
Vacuum pump SV 40	34063

Order no.	
Foam sensor with distribution head cpl.	40507



### Splash protector

A polycarbonate splash protector mounted on the edge of the bath, protects against splashing from the rotating flask.



### Vapor duct with frit P3

This special vapor duct protects the condenser assembly against powder and dust during the drying process.



### Flask handler crane

With the folding crane and the 50 L manual flask handler (included) the glass flasks are conveniently and safely secured and transported.



### Evaporating and drying flasks

With large opening (inner diameter 120 mm) for easy sample taking and cleaning.

Order no.	
Splash protector for R-250	41420

Order no.	
Vapor duct with frit P3	41100

Order no.	
Flask handler crane for R-250	41494
Flask handler crane for R-250 EX	41493

Order no.	
20-liter flask	41432
50-liter flask	41339
20-liter drying flask	41393
50-liter drying flask	41394



### Manual flask handler

Assists in easy mounting and removal of large volume flask along with safe transport and sturdy base for positioning.



### Vacuum pump MD4C EX

ATEX compliant diaphragm vacuum pump, oil-free and without sources of ignition for the EX area. The flow rate is 3.7 m<sup>3</sup>/h with an ultimate pressure of 3 mbar.

Order no.	
Manual flask handler for 20-liter flasks on R-250	41410
Manual flask handler for 50-liter flasks on R-250	41414

Order no.	
Vacuum pump MD4C EX	40197



**Only accessories with this symbol are to be used in EX zones.**

For optimized distillations, the Rotavapor is paired with supporting equipment for efficient cooling, vacuum control and gradients.



In addition to vacuum pumps, controllers and sensors, Buchi offers a compatible R-220 SE recirculating chiller. With a cooling performance of 2500 W at 15 °C it's an excellent energy efficient alternative to cooling water.

Dimensions (WxHxD): 45 x 110 x 55 cm

Weight: 72 kg

Tank volume: 20 liters

Control range: -10 °C ... 40 °C

Power supply: 230 V/50 Hz

Order number

40300

## Customized solutions

If specific applications or physical restrictions do not permit the usage of our standard instruments, Buchi can accommodate a wide variety of custom designs; ranging from a special lid for the evaporating flask to a completely modified rotary evaporator.

### Example 1:

#### R-250 with low height

An R-250 with reflux required a reduced height of 190 cm to fit into a laboratory with spacial limitations.

By modifying individual glass components Buchi was able to meet the customer's special requirement.

### Example 2:

#### Lid for sealing the evaporating flask

A custom polypropylene lid with EPDM O-ring for sealing the evaporating flask is used to protect the product from contaminants.

Order number

42895

### Example 3:


#### Heating bath with suspended cooling coil

The heating bath was equipped with a suspended cooling to expedite cool down times and eliminate the need to change the bath's contents. The efficient design has two connections for a recirculating chiller. This was used for recrystallization.




# Order information



## The new Rotavapor® R-220 SE



R-220 SE for continuous distillation	Order no.
Version 230V 3.6kW Heating bath oil and water	118114C291
Version 200V 3.6kW Heating bath oil and water	118124C291
Version 400V 3.6kW Heating bath oil and water	118164C291
Version 400V 4.2 kW Water bath	118194C291





## Rotavapor® R-220 EX





## Rotavapor® R-250 EX



### R-250 EX for continuous distillation

	Order no.
Version T3 up to 150 °C	8610D580
Version T4 up to 135 °C	8620D580

## Application specific suggestions

Standard setup R-220 SE		Order no.
R-220 SE with:		
	Vacuum pump V-710	72000
	Manual flask handler	41400
	Cooling water valve	41191
	Trolley	41257
Additional for performance optimization		Order no.
	Coolant sensor	11055988
	Level sensors	11056082
For maximum safety		Order no.
	Flow monitor	11055971
	Coolant sensor	11055988
	Safety shield cpl.	Depending on glass assy.
For GMP		Order no.
	USB module	11056175
	Coolant sensor	11055988

# Technical data

**NEW**

	R-220 SE	R-250
<b>Distillation capacity (ethanol)</b>	Up to 13.5 liters/hour 3600 W heater Up to 16 liters/hour 4200 W heater	Up to 19.5 liters/hour
<b>Electrical</b>		
Voltage	200V, 220–240V (1P, N, G) / 400V (3P, N, G)	400–440 V (3P, N, E)
Frequency	50–60 Hz	50–60 Hz
Power consumption	4200 W	7500 W–9100 W
<b>Dimensions</b>		
Height	1430 mm–1750 mm, depending on glass assy.	2100 mm–2300 mm, depending on glass assy.
Width	1200 mm	1410 mm
Depth	510 mm	780 mm
Bath diameter (inner)	430 mm	560 mm
<b>Weight</b>	Approx. 65 kg without glass	Approx. 160 kg without glass
<b>Approvals</b>		
CE lable	Yes	Yes
EX lable	–	–
<b>Safety</b>		
Bath over temperature protection	– Separate electronic monitoring circuit, mechanical reset – Error if the actual temperature is 15 °C above the setpoint	
Bath level monitoring	Not installed	Not installed
Rotation	Soft start	Soft start
Behavior in case of error	– Bath lowered, heater off, rotation off – Indication of the error – Reset possible	– Bath lowered, heater off, rotation off – Indication of the error – Reset possible
<b>Display</b>		
Bath temperature	Yes (1 °C steps)	Yes (1 °C steps)
Vapor temperature	Yes (1 °C steps)	Yes (1 °C steps)
Setting for the rotation speed	Yes (1 rpm steps)	Yes (1 rpm steps)
Setting for the bath temperature	Yes (1 °C steps)	Yes (1 °C steps)
Actual vacuum	Yes (1 mbar steps)	Optional vacuum controller
Setting for the vacuum	Yes (1 mbar steps)	Optional vacuum controller
<b>Cooling</b>		
Cooling water flow rate	120–200 liters/hour (variable due to integ. valve)	200–400 liters/hour (variable due to integ. valve)
Limitations	Max. 2.7 bar abs. pulse-free	Max. 2.7 bar abs. pulse-free
<b>Heating</b>		
Medium	Water or oil	Water or oil
Temperature range	Room temperature up to 180°	Room temperature up to 180°
Heating power	230 V 1-phase; 3600 W (3 W/cm <sup>2</sup> ) optional 4200 W (only for water)	400 V 3-phase; 6600 W (3 W/cm <sup>2</sup> )
Regulation precision	+/- 2 °C	+/- 2 °C
<b>Rotation</b>		
Motor	200 V 1-phase; 0.6 A at 50 Hz; 1Nm	400 V 3-phase; 2 A at 50 Hz; 3 Nm
Control	Electronic	Frequency inverter
Rotation speed	5–140 rpm	5–100 rpm
Accuracy	+/- 1 upm at 5 rpm up to +/- 5 rpm at 150 rpm	+/- 1 upm at 5 rpm up to +/- 5 rpm at 100 rpm
<b>Materials</b>		
Housing	Stainless steel 1.4301	Stainless steel 1.4301
Gearbox head	Aluminum (3.2373)	Cast aluminum (3.2373)
Coating	Powder-coated with Epoxit (EPX)	Powder-coated with Epoxit (EPX)
Bath construction	Stainless steel 1.4404	Stainless steel 1.4404
Heating element	Stainless steel 1.4404	Stainless steel 1.4404
Glass	Borosilicate 3.3	Borosilicate 3.3
Product contacting parts	FDA compliant materials	FDA compliant materials
<b>Pressure in the glass assembly</b>	0 mbar up to ambient pressure	0 mbar up to ambient pressure
Tightness	< 1 mbar/min.	< 1 mbar/min.
<b>Sensors</b>		
Vapor temperature	PT-1000, 2-wire	PT-1000, 2-wire
Bath temperature	PT-1000, 2-wire	PT-1000, 2-wire
Bath level sensor	Not installed	Not installed
Pressure	Ceramic, capacitive	Optional vacuum controller

	R-220 EX	R-250 EX
<b>Distillation capacity</b>	Up to 13.5 liters/hour	Up to 19.5 liters/hour
<b>Electrical</b>		
Voltage	200 V, 220–240 V (1P, N, G) / 400 V (3P, N, G)	400–440 V (3P, N, E)
Frequency	50–60 Hz	50–60 Hz
Power consumption	4200 W	7500 W–9100 W
<b>Dimensions</b>		
Height	1430 mm–1750 mm, depending on glass assy.	2100 mm – 2300 mm, depending on glass assy.
Width	1200 mm	1410 mm
Depth	660 mm	850 mm
Bath diameter (inner)	430 mm	560 mm
<b>Weight</b>	Approx. 75 kg without glass	Approx. 250 kg without glass
<b>Approvals</b>		
CE lable	Yes	Yes
EX lable for Europe	EEx de ia [ia] IIC T3 and T4	EEx de ia [ia] IIC T3 and T4
<b>Safety</b>		
Bath over temperature protection	– Separate electronic monitoring circuit, mechanical reset – Error if the actual temperature is 15 °C above the setpoint	
Bath level monitoring	Prevents dry operation of the bath	Prevents dry operation of the bath
Rotation	Soft start	Soft start
Behavior in case of error	– Heater off – Indication of the error – Reset possible	– Bath lowered, heater off, rotation off – Indication of the error – Reset possible
<b>Display</b>		
Bath temperature	Yes (1 °C steps)	Yes (1 °C steps)
Vapor temperature	Yes (1 °C steps)	Yes (1 °C steps)
Setting for the rotation speed	Yes (no display)	Yes (1 rpm steps)
Setting for the bath temperature	Yes (1 °C steps)	Yes (1 °C steps)
Actual vacuum	Yes (1 mbar steps)	Yes (1 mbar steps)
Setting for the vacuum	Yes (1 mbar steps)	Yes (1 mbar steps)
<b>Cooling</b>		
Cooling water flow rate	120–200 liters/hour (variable due to integrated valve)	200–400 liters/hour (variable due to integrated valve)
Limitations	Max. 2.7 bar abs. pulse-free	Max. 2.7 bar abs. pulse-free
<b>Heating</b>		
Medium	Water or oil	Water or oil
Temperature range	Room temperature up to 135 °C (T4) or 150 °C (T3)	Room temperature up to 135 °C (T4) or 150 °C (T3)
Heating power	230 V 1-phase; 3600 W (3 W/cm <sup>2</sup> )	400 V 3-phase; 6600 W (3 W/cm <sup>2</sup> )
Regulation precision	+/- 2 °C	+/- 2 °C
<b>Rotation</b>		
Motor	Pneumatic	400 V 3-phase; 2 A to 50 Hz; 3 Nm
Control	Continuous with precision regulator valve	Frequency inverter
Rotation speed	0–100 rpm	5–100 rpm
Accuracy	No display	+/- 1 rpm at 5 rpm up to +/- 5 rpm at 100 rpm
<b>Materials</b>		
Housing	Stainless steel 1.4301	Stainless steel 1.4301
Gearbox head	Aluminum (3.2373)	Cast aluminum (3.2373)
Coating	Powder-coated with Epoxit (EPX)	Powder-coated with Epoxit (EPX)
Bath construction	Stainless steel 1.4404	Stainless steel 1.4404
Heating element	Stainless steel 1.4404	Stainless steel 1.4404
Glass	Borosilicate 3.3	Borosilicate 3.3
Product contacting parts	FDA-compliant materials	FDA compliant materials
<b>Pressure in the glass assembly</b>	0 mbar up to ambient pressure	0 mbar up to ambient pressure
Tightness	< 1 mbar/min.	< 1 mbar/min.
<b>Sensors</b>		
Vapor temperature	PT-1000, 2-wire	PT-1000, 2-wire
Bath temperature	PT-1000, 2-wire	PT-1000, 2-wire
Bath level sensor	Ultrasonic	Ultrasonic
Pressure	Piezoresistive measuring cell; Stainless steel 1.4435	Piezoresistive measuring cell; Stainless steel 1.4435

# Buchi products and services

**Buchi provides laboratory equipment for research and development as well as analytical instruments for quality control – worldwide. All our products embody our philosophy of “Quality in your hands”. We meet the needs of our customers to the highest degree by developing products and solutions that are robust, well designed, and easy to use.**



## Buchi products and services:



BÜCHI Labortechnik AG  
Postfach  
9230 Flawil 1  
Switzerland  
T +41 71 394 63 63  
F +41 71 394 65 65  
buchi@buchi.com  
www.buchi.com

BÜCHI Labortechnik GmbH  
Postfach 10 03 51  
45003 Essen  
Germany  
Freecall 0800 414 0 414  
T +49 201 747 490  
F +49 201 237 082  
deutschland@buchi.com  
www.buechigmbh.de

BUCHI Sarl  
5, rue du Pont des Halles  
Z.A. du Delta  
94656 Rungis Cedex  
France  
T +33 1 56 70 62 50  
F +33 1 46 86 00 31  
france@buchi.com  
www.buchi.fr

BÜCHI Italia s.r.l.  
Centro Direzionale, Milano Fiori  
Pal. A-4, Strada 4  
20090 Assago (MI)  
Italia  
T +39 02 824 50 11  
F +39 02 57 51 28 55  
italia@buchi.com  
www.buchi.it

BÜCHI Labortechnik GmbH  
Branch Office Netherlands  
Postbus 142  
3340 AC Hendrik-Ido-Ambacht  
The Netherlands  
T +31 78 684 94 29  
F +31 78 684 94 30  
netherlands@buchi.com  
www.buchi.nl

BUCHI Corporation  
19 Lukens Drive, Suite 400  
New Castle  
Delaware 19720  
USA  
T +1 302 652 3000  
F +1 302 652 8777  
Toll Free: +1 877 692 8244  
us-sales@buchi.com  
www.mybuchi.com

Nihon BUCHI K.K.  
3F IMON Bldg.,  
2-7-17 Ikenohata, Taito-ku,  
Tokyo 110-0008  
Japan  
T +81 3 3821 4777  
F +81 3 3821 4555  
nihon@buchi.com  
www.nihon-buchi.jp

BUCHI India  
Private Ltd.  
201, Magnum Opus  
Shantinagar Industrial Area  
Vakola, Santacruz (East)  
Mumbai 400 055,  
India  
T +91 22 667 18983 / 84 / 85  
F +91 22 667 18986  
www.buchi.com

BUCHI UK Ltd  
5 Whitegate Business Centre  
Jardine Way  
Chadderton  
Oldham OL9 9QL  
United Kingdom  
T +44 161 633 1000  
F +44 161 633 1007  
uk@buchi.com  
www.buchi.co.uk

BUCHI Hong Kong Ltd.  
14<sup>th</sup> Floor South China Building  
1-3 Wyndham Street  
Central, Hong Kong  
China  
T +852 2389 2772  
F +852 2389 2774  
china@buchi.com  
www.buchi.com.cn

BUCHI Shanghai  
570 West Huaihai Road  
NO.38 Chuangyi Street  
500052 Shanghai  
China  
T +86 21 6280 3366  
F +86 21 5230 8821  
china@buchi.com  
www.buchi.com.cn

BUCHI (Thailand) Ltd.,  
77/121, Sin Sathon Tower,  
28th FL, Unit C  
Krunghonburi Rd.  
Klongtongsai, Klongsan  
Bangkok 10600  
Thailand  
T +66 2 862 08 51  
F +66 2 862 08 54  
bacc@buchi.com  
www.buchi.com

We are represented by more than 100 distribution partners worldwide. Find your local representative at [www.buchi.com](http://www.buchi.com)

**Quality in your hands**

