

## **Non-Stirred Pressure Vessels**

### **General Purpose Non-stirred Reactors and Pressure Vessels**

Parr offers a broad selection of non-stirred pressure vessels in convenient styles, sizes and pressure ranges for many laboratory uses. Recent additions to this line have included new sizes, new temperature and pressure combinations and new self-sealing closures, all designed to provide convenient vessels for laboratory procedures which must be performed under pressure at elevated temperatures. For more information, please read our 'Selection Guide to Parr Pressure Vessels.'



#### **Series 4600 - 4620 GP 4.00 Inch Diameter**

Sizes: 1000mL, 2000mL



#### **Series 4632 Pressure Extraction Vessel 4.00 Inch Diameter**

Size: 2000 mL



#### **Series 4650 2.50 Inch Inside Diameter HP/HT**

Sizes: 250mL, 500mL, 1000mL



#### **Series 4660 6.00 Inch Inside Diameter GP**

Sizes: 3.75L, 7.5L



**Series 4670 5.50 Inch Inside Diameter HP/HT**

Sizes: 3750mL, 5800mL



**Series 4680 3.75 Inch Inside Diameter HP/HT**

Sizes: 1000mL, 1800mL



**Series 4676 9.50 Inch Diameter GP**

Sizes: 10000, 20000mL



**Series 4750 1.50 Inch Inside diameter GP**

Sizes: 125mL, 200mL



**Series 4760 - 4774 2.50 Inch Inside Diameter GP**

Sizes: 160mL, 600mL



**Series 4790 Micro GP**

Sizes: 25 mL, 50mL, 100 mL



**Series 4703 - 4714 Screw Cap Vessels 1.00 Inch Inside Diameter**

Sizes: 22mL, 45mL



**Series 4740 1.00 Inch Inside Diameter HP/HT**

Sizes: 21mL, 71mL